

Product Datasheet

Adoprazine (orb1710486)

Description	Adoprazine
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{24}H_{24}FN_3O_2$
Purity	99.44%
MW	405.46
CAS Number	222551-17-9
SMILES	<chem>C(N1CCN(C2=C3C(=CC=C2)OCCO3)CC1)C=4C=C(C=NC4)C5=CC=C(F)C=C5</chem>
Expiration Date	12 months from date of receipt.