

Product Datasheet

A2AR-agonist-1 (orb1710514)

Description	A2AR-agonist-1
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{20}H_{22}N_6O_4$
Purity	99.82%
MW	410.43
CAS Number	41552-95-8
SMILES	<chem>O[C@H]1[C@H](N2C=3C(N=C2)=C(NCCC=4C=5C(NC4)=CC=CC5)N=C N3)O[C@H](CO)[C@H]1O</chem>
Expiration Date	12 months from date of receipt.