

Product Datasheet

Recombinant Tau(Phospho-Thr217) Rabbit mAb(G122) Antibody (orb1710795)

Description	Recombinant Tau(Phospho-Thr217) Rabbit mAb(G122)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit
Uniprot ID	P10636
Expiration Date	12 months from date of receipt.