

Product Datasheet

Canine PDGF R beta / CD140b Protein, His Tag (orb1710808)

Description	Canine PDGF R beta Protein, His Tag is expressed from human 293 cells (HEK293). It contains AA Leu 3...
Reactivity	Canine
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PDGF R beta
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 58.0 kDa. The protein migrates as 70-105 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Canine PDGF R beta Protein, His Tag (orb1710808) is expressed from human 293 cells (HEK293). It contains AA Leu 33 - Val 532 (Accession # Q6QNF3-1).
MW	58.0 kDa
Uniprot ID	Q6QNF3
Expiration Date	6 months from date of receipt.