



Product Datasheet

Interferon alpha/beta receptor 1 Antibody (orb1710816)


Description	Interferon alpha/beta receptor 1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP, WB
Immunogen	Residues 27-436 ECD)
Target	Interferon alpha/beta receptor 1
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	400 µg/ml
Storage	Shelf life: 1 year from date of receipt. Storage: 2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	P17181
Dilution Range	WB: 1:1,000, IP: 15 µl/1 mg lysate, IHC: 1:1000 to 1:4000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human Interferon alpha/beta...



Detection of human Interferon alpha/beta...



Detection of human Interferon alpha/beta...