



## Product Datasheet

### GLI1 Antibody (orb1710825)

<b>Description</b>	GLI1 Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IP, WB
<b>Immunogen</b>	Between 1056 and C-terminus
<b>Target</b>	GLI1
<b>Preservatives</b>	Borate Buffered Saline (BBS) pH 8.2 with 0.1% rAlbumin and 0.09% Sodium Azide
<b>Form/Appearance</b>	Whole IgG
<b>Concentration</b>	100 µg/ml
<b>Storage</b>	Shelf life: 1 year from date of receipt. Storage: 2 - 8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Format: Whole IgG
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<b>P08151</b>
<b>Dilution Range</b>	WB: 1:1,000, IP: 20 µl/1 mg lysate
<b>Expiration Date</b>	12 months from date of receipt.



Detection of human GLI1 by western blot....



Detection of human GLI1 by western blot ...