

Product Datasheet

TAG-72 / CA72.4 (Tumor-Associated Glycoprotein); Clone SPM536 (Concentrate) Antibody (orb1712088)

Description	TAG-72 / CA72.4 (Tumor-Associated Glycoprotein); Clone SPM536 (Concentrate) Antibody
Species/Host	Mouse
Reactivity	Bovine, Canine, Human, Rat
Conjugation	Unconjugated
Immunogen	Purified TAG-72 protein.
Note	For research use only
Application notes	Jurkat cells. Breast or lung carcinoma.
Isotype	IgG1, kappa.
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.