

Product Datasheet

ER beta-1 (Estrogen Receptor beta-1) Antibody (orb1712376)

Catalog Number	orb1712376
Category	Antibodies
Description	ER beta-1 (Estrogen Receptor beta-1); Clone ERb455 (Concentrate)
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG2a
Conjugation	Unconjugated
Reactivity	Equine, Human, Monkey, Mouse, Porcine, Rat, Sheep
Immunogen	C-terminus fragment of recombinant human estrogen receptor beta protein
Application notes	Ovarian, breast or prostate carcinoma.
Clone Number	ERb455
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)