

Product Datasheet

CD1a / HTA1 (Mature Langerhans Cells Marker); Clone O10 (Concentrate) Antibody (orb1712609)

Description	CD1a / HTA1 (Mature Langerhans Cells Marker); Clone O10 (Concentrate) Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Immunogen	Human thymus cells
Note	For research use only
Application notes	MOLT-4 cells. Paracortex in a tonsil or a reactive lymph node.
Isotype	IgG1, kappa
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.