

## Product Datasheet

# CD1a / HTA1 (Mature Langerhans Cells Marker); Clone O10 & C1A/711 (Concentrate) Antibody (orb1712610)

<b>Description</b>	CD1a / HTA1 (Mature Langerhans Cells Marker); Clone O10 & C1A/711 (Concentrate) Antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Immunogen</b>	Human thymus cells (O10); Recombinant human CD1a protein (C1A/711)
<b>Note</b>	For research use only
<b>Application notes</b>	MOLT-4 cells. Paracortex in a tonsil or a reactive lymph node.
<b>Isotype</b>	IgG1, kappa (O10) & IgG1, kappa (C1A/711)
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.