

Product Datasheet

TAG-72 / CA72.4 (Tumor-Associated Glycoprotein); Clone SPM536 (Concentrate) Antibody (orb1712769)

| | |
|--------------------------|--|
| Description | TAG-72 / CA72.4 (Tumor-Associated Glycoprotein); Clone SPM536 (Concentrate) Antibody |
| Species/Host | Mouse |
| Reactivity | Bovine, Canine, Human, Rat |
| Conjugation | Unconjugated |
| Immunogen | Purified TAG-72 protein. |
| Note | For research use only |
| Application notes | Jurkat cells. Breast or lung carcinoma. |
| Isotype | IgG1, kappa. |
| Clonality | Monoclonal |
| Expiration Date | 12 months from date of receipt. |