

Product Datasheet

TAG-72 / CA72.4 (Tumor-Associated Glycoprotein); Clone SPM148 (Concentrate) Antibody (orb1712770)

Description	TAG-72 / CA72.4 (Tumor-Associated Glycoprotein); Clone SPM148 (Concentrate) Antibody
Species/Host	Mouse
Reactivity	Bovine, Canine, Hamster, Human, Rat
Conjugation	Unconjugated
Immunogen	Membrane-enriched fraction of a human breast carcinoma liver metastasis.
Note	For research use only
Application notes	Jurkat cells. Breast or lung carcinoma.
Isotype	IgG1, kappa.
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.