

Product Datasheet

TAG-72 (Tumor-Associated Glycoprotein); Clone CA72/733 (Concentrate) Antibody (orb1712772)

Description	TAG-72 (Tumor-Associated Glycoprotein); Clone CA72/733 (Concentrate) Antibody
Species/Host	Mouse
Reactivity	Bovine, Canine, Human, Rat
Conjugation	Unconjugated
Immunogen	Recombinant human TAG-72 protein
Note	For research use only
Application notes	Jurkat cells, breast or lung carcinoma.
Isotype	IgG1, kappa
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.