

Product Datasheet

ER beta-1 (Estrogen Receptor beta-1) Antibody (orb1713064)

Description	ER beta-1 (Estrogen Receptor beta-1); Clone ERb455 (Concentrate)
Species/Host	Mouse
Reactivity	Equine, Human, Monkey, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Immunogen	C-terminus fragment of recombinant human estrogen receptor beta protein
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Ovarian, breast or prostate carcinoma.
Isotype	IgG2a
Clonality	Monoclonal
Clone Number	ERb455
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)