

Product Datasheet

CD1a / HTA1 (Mature Langerhans Cells Marker); Clone O10 & C1A/711 (Concentrate) Antibody (orb1713306)

Description	CD1a / HTA1 (Mature Langerhans Cells Marker); Clone O10 & C1A/711 (Concentrate) Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Immunogen	Human thymus cells (O10); Recombinant human CD1a protein (C1A/711)
Note	For research use only
Application notes	MOLT-4 cells. Paracortex in a tonsil or a reactive lymph node.
Isotype	IgG1, kappa (O10) & IgG1, kappa (C1A/711)
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.