

Product Datasheet

Bcl-X (Apoptosis Marker); Clone BX006 & 2H12 (Concentrate) Antibody (orb1713366)

Description	Bcl-X (Apoptosis Marker); Clone BX006 & 2H12 (Concentrate) Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Immunogen	A synthetic peptide, aa 3-14 (Cys-QSNRELVVDFLS) of human Bcl-X protein
Note	For research use only
Application notes	Jurkat, K562, HL-60, or HeLa Cells. Reed-Sternberg cells in Hodgkin's lymphoma.
Isotype	IgG2a (BX006) & IgG2a (2H12)
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.