

## Product Datasheet

# Bcl-X (Apoptosis Marker); Clone 2H12 (Concentrate) Antibody (orb1713367)

<b>Description</b>	Bcl-X (Apoptosis Marker); Clone 2H12 (Concentrate) Antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Porcine, Rat
<b>Conjugation</b>	Unconjugated
<b>Immunogen</b>	A synthetic peptide, aa 3-14 (Cys-QSNRELVVDFLS) of human Bcl-X protein
<b>Note</b>	For research use only
<b>Application notes</b>	Jurkat, K562, HL-60, or HeLa Cells. Reed-Sternberg cells in Hodgkin's lymphoma.
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.