

Product Datasheet

Bexarotene (orb1713587)

Description	Retinoid RXR agonist
Form/Appearance	White powder
Storage	-20°C
Note	For research use only
Formula	$C_{24}H_{28}O_2$
Purity	> 98% (TLC); NMR (Conforms)
Source	Synthetic
MW	348.5
CAS Number	153559-49-0
Solubility (25°C)	Soluble in DMSO (60 mg/ml) or ethanol (warm, 10 mg/ml)
Hazard Information	Classification: Not a hazardous substance or mixture. Safety Phrases: S22 - Do not breathe dust. S24/25 - Avoid contact with skin and eyes. S36/37/39 - Wear suitable protective clothing, gloves and eye/face protection.
Expiration Date	12 months from date of receipt.