

Product Datasheet

Lestaurtinib (orb1713728)

Description	FLT3, TrkA kinase inhibitor
Form/Appearance	White powder
Storage	-20°C
Note	For research use only
Formula	$C_{26}H_{21}N_3O_4$
Purity	≥98% (HPLC)
Source	Synthetic
MW	439.5
CAS Number	111358-88-4
Solubility (25°C)	Soluble to 100 mM in DMSO and to 25 mM in ethanol
Hazard Information	Classification: Not a hazardous substance or mixture. Safety Phrases: S22 - Do not breathe dust. S24/25 - Avoid contact with skin and eyes. S36/37/39 - Wear suitable protective clothing, gloves and eye/face protection.
Expiration Date	12 months from date of receipt.