

Product Datasheet

BzATP (orb1713804)

Description	Prototypic P2X7 receptor agonist
Form/Appearance	Off white solid
Storage	-20°C
Note	For research use only
Formula	$C_{24}H_{24}N_5O_{15}P_3 \cdot C_{18}H_{45}N_3$
Purity	> 97.7%
Source	Synthetic
MW	1018.97
CAS Number	112898-15-4
Solubility (25°C)	Soluble in water.
Hazard Information	Classification: Not WHMIS controlled. Safety Phrases: S22 - Do not breathe dust. S24/25 - Avoid contact with skin and eyes. S36/37/39 - Wear suitable protective clothing, gloves and eye/face protection.
Expiration Date	12 months from date of receipt.