

Product Datasheet

HSP60 Protein (orb1714300)

Description	Active Human Recombinant HSP60 Protein
Reactivity	Human
Tested Applications	ATPase Activity Assay, ELISA, FA, SDS-PAGE, WB
Target	HSP60
Preservatives	20mM Phosphate Buffer, 150mM NaCl, 10% glycerol
Concentration	Lot/batch specific. See included datasheet.
Storage	-20°C
Tag	His tag
Note	For research use only
Application notes	This product has been certified > 90% pure using SDS-PAGE analysis. The protein has ATPase activity at the time of manufacture of 3.6µM phosphate liberated/hr/µg protein in a 200µl reaction at 37°C (pH7.5) in the presence of 20ul of 1mM ATP using a Malachite Green assay.
Purity	> 90%
Source	Recombinant
MW	~60 kDa
Uniprot ID	P10809
Expiration Date	6 months from date of receipt.