
Product Datasheet

HSP90 beta Protein (orb1714303)

Description

Active Human Recombinant HSP90 beta Full Length Protein

Reactivity

Human

Tested Applications

CoIP, ELISA, FA, SBR, SDS-PAGE, WB

Target

HSP90 beta

Preservatives

20mM Tris, pH 7.5, 175 mM NaCl, 0.1 mM EDTA, 10% glycerol, 1 mM DTT

Concentration

Lot/batch specific. See included datasheet.

Storage

-20°C

Tag

No tag

Note

For research use only

Application notes

This product has been certified > 90% pure using SDS-PAGE analysis.

Protein Sequence

MPEEVHHGEEEVETFAFQAEIAQLMSLIINTFYNSKEIFLRELISNASDALDKIRYES
 LTDPSKLDGKELKIDIIPNPQERTLTLVDTGIGMTKADLNNLGTIAKSGTKAFME
 ALQAGADISMIGQFGVGFYSAYLVAEKVVITKHNDDEQYAWESSAGGSFTVRA
 DHGEPGRGTVILHLKEDQTEYLEERRVKEVVKKHSQFIGYPITLYLEKEREKEIS
 DDEAEKEEKEEEDKDEEKPKIEDVGSDEEDDSGKDKKKKTKKIKEKYIDQEE
 LNKTPIWTRNPDDITQEEYGEFYKSLTNDWEDHLAVKHFSVEGQLEFRALLFIP
 RRAPFDLFENKKKNNIKLYVRRVFMDSCELIPEYLNFIKRVVDSDELPLNISRE
 MLQQSKILKVIKKNIVKKCLELSELAEDKENYKFFYAFSKNLKLGIHEDSTNRRR
 LSELLRYHTSQSGDEMTSLSEYVSRMKETQKSIYYITGESKEQVANSFVERVRK
 RGFVVMTEPIDEYCVQLKEFDGKSLVSVTKEGLELPEDEEKKKMEESKAKF
 ENLCKLMKEILDKKVEKVTISNRLVSSPCCIVTSTYGWTANMERIMKAQALRDN
 TMGYMMAKKHLEINPDHPIVETLRQKAEADKNDKAVKDLVLLFETALLSSGFSL
 EDPQTHSNRIYRMIKGLGLIDEDEVAEEPNAAVPDEIPPLEGDEDEDASRMEEVD

Purity

> 85%

Source

Recombinant

MW

~90 kDa

Uniprot ID
P08238
Expiration Date

6 months from date of receipt.