

Product Datasheet

anti- A29L antibody (orb1714623)

Description	anti- A29L antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Monkeypox Virus (MPXV) Protein A29L
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	Q77HM6
Dilution Range	WB: 1:1000-1:5000; ELISA: 1:5000 - 1:10000
Expiration Date	12 months from date of receipt.