


## Product Datasheet

### T3(Triiodothyronine) ELISA Kit (orb1714641)

<b>Description</b>	T3(Triiodothyronine) ELISA Kit
<b>Reactivity</b>	All
<b>Range</b>	0.078-5ng/ml
<b>Target</b>	T3
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Sensitivity</b>	0.047ng/ml
<b>Expiration Date</b>	Please enquire.



T3(Triiodothyronine)  
ELISA Kit