


## Product Datasheet

# Rabbit FSH(Follicle-Stimulating Hormone)ELISA Kit (orb1714648)

<b>Description</b>	Rabbit FSH(Follicle-Stimulating Hormone)ELISA Kit
<b>Reactivity</b>	Rabbit
<b>Range</b>	1.563-100mIU/ml
<b>Target</b>	FSH
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<a href="#">Q6IY73</a>
<b>Sensitivity</b>	0.938mIU/ml
<b>Expiration Date</b>	Please enquire.



Rabbit FSH(Follicle-Stimulating Hormone)...