


Product Datasheet

Rat HS-GH(High sensitive Growth Hormone) Accquant ELISA Kit (orb1714649)

Description	Rat HS-GH(High sensitive Growth Hormone) Accquant ELISA Kit
Reactivity	Rat
Range	3.125-200pg/ml
Target	HS-GH
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P01244
Sensitivity	1.875pg/ml
Expiration Date	Please enquire.



Rat HS-GH(High sensitive Growth Hormone)...