


Product Datasheet

Porcine C5a (Complement Component 5a) ELISA Kit (orb1714664)

| | |
|------------------------|--|
| Description | Porcine C5a (Complement Component 5a) ELISA Kit |
| Reactivity | Porcine |
| Range | 78.125-5000pg/ml |
| Target | C5a |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Sample Types | Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples |
| Uniprot ID | P01032 |
| Sensitivity | 46.875pg/ml |
| Expiration Date | Please enquire. |



Porcine C5a (Complement Component 5a) EL...