


Product Datasheet

Human PSMD8(26S proteasome non-ATPase regulatory subunit 8)ELISA Kit (orb1714719)

Description	Human PSMD8(26S proteasome non-ATPase regulatory subunit 8)ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	PSMD8
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P48556
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
PSMD8(26S
proteasome
non-ATPase
re...