


Product Datasheet

Human IL36RN(Interleukin-36 receptor antagonist protein)ELISA Kit (orb1714721)

Description	Human IL36RN(Interleukin-36 receptor antagonist protein)ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	IL36RN
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9UBH0
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
IL36RN(Interleukin-
36 receptor ant...