


Product Datasheet

Human TCC C5b-9(Terminal Complement Complex C5b-9)ELISA Kit (orb1714729)

Description	Human TCC C5b-9(Terminal Complement Complex C5b-9)ELISA Kit
Reactivity	Human
Range	12.5-800ng/ml
Target	TCC C5b-9
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	7.5ng/ml
Expiration Date	Please enquire.



Human TCC C5b-9(Terminal Complement Comp...