


Product Datasheet

Human GLP2(Glucagon Like Peptide 2) ELISA Kit (orb1714734)

Description	Human GLP2(Glucagon Like Peptide 2) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	GLP2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P01275
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
GLP2(Glucagon
Like Peptide 2)
ELIS...