


Product Datasheet

Rat MKRN-1(Macorin Ring Zinc Finger Protein-1)ELISA kit (orb1715030)

Description	Rat MKRN-1(Macorin Ring Zinc Finger Protein-1)ELISA kit
Reactivity	Rat
Range	31.25-2000pg/ml
Target	MKRN-1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	A0A8I6ALB9
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat MKRN-1(Macorin Ring Zinc Finger Prot...