


Product Datasheet

Human PSMD2(26S proteasome non-ATPase regulatory subunit 2)ELISA Kit (orb1715058)

Description	Human PSMD2(26S proteasome non-ATPase regulatory subunit 2)ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	PSMD2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q13200
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human PSMD2(26S proteasome non-ATPase re...