


## Product Datasheet

# Human cpeb4(Cytoplasmic polyadenylation element-binding protein 4)ELISA Kit (orb1715059)

<b>Description</b>	Human cpeb4(Cytoplasmic polyadenylation element-binding protein 4)ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	cpeb4
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q17RY0</b>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
cpeb4(Cytoplasmic  
polyadenylation ...