


## Product Datasheet

### Human PSMA7(Proteasome subunit alpha type-7)ELISA Kit (orb1715078)

<b>Description</b>	Human PSMA7(Proteasome subunit alpha type-7)ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.25-2000pg/ml
<b>Target</b>	PSMA7
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>O14818</b>
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
PSMA7(Proteasome  
subunit alpha typ...