


## Product Datasheet

# Human anti ET-1(anti-Endothelin 1 IgG Antibody)ELISA Kit (orb1715085)

<b>Description</b>	Human anti ET-1(anti-Endothelin 1 IgG Antibody)ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	3.125-200ng/ml
<b>Target</b>	anti ET-1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sensitivity</b>	1.875ng/ml
<b>Expiration Date</b>	Please enquire.



Human anti ET-1(anti-Endothelin 1 IgG An...