


Product Datasheet

Rabbit GHRL(Ghrelin)ELISA Kit (orb1715105)

Description	Rabbit GHRL(Ghrelin)ELISA Kit
Reactivity	Rabbit
Range	3.125-200pg/ml
Target	Ghrelin
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	A5JSY6
Sensitivity	1.875pg/ml
Expiration Date	Please enquire.



Rabbit
GHRL(Ghrelin)
ELISA Kit