


## Product Datasheet

### Goat Gc(Glucagon)ELISA Kit (orb1715118)

<b>Description</b>	Goat Gc(Glucagon)ELISA Kit
<b>Reactivity</b>	Goat
<b>Range</b>	31.25-2000pg/ml
<b>Target</b>	Glucagon
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">A0A452FHY7</a>
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.



Goat  
Gc(Glucagon)  
ELISA Kit