


Product Datasheet

Rat TSH(Thyroid Stimulating Hormone)ELISA Kit (orb1715124)

Description	Rat TSH(Thyroid Stimulating Hormone)ELISA Kit
Reactivity	Rat
Range	0.625-40ng/ml
Target	TSH
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P04652
Sensitivity	0.375ng/ml
Expiration Date	Please enquire.



Rat
TSH(Thyroid
Stimulating
Hormone)
ELI...