


## Product Datasheet

# Human Anti-KLHL12(anti-Kelch-like protein 12 antibody)ELISA Kit (orb1715125)

<b>Description</b>	Human Anti-KLHL12(anti-Kelch-like protein 12 antibody)ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.563-100ng/ml
<b>Target</b>	Anti-KLHL12
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q53G59</a>
<b>Sensitivity</b>	0.938ng/ml
<b>Expiration Date</b>	Please enquire.



Human Anti-KLHL12(anti-Kelch-like protei...