


## Product Datasheet

# Human PSMB4(Proteasome subunit beta type-4) ELISA Kit (orb1715134)

<b>Description</b>	Human PSMB4(Proteasome subunit beta type-4) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/ml
<b>Target</b>	PSMB4
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P28070</b>
<b>Sensitivity</b>	0.094ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
PSMB4(Proteasome  
subunit beta  
type...