


## Product Datasheet

# Human YTHDF3(YTH domain-containing family protein 3) ELISA Kit (orb1715167)

<b>Description</b>	Human YTHDF3(YTH domain-containing family protein 3) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.25-2000pg/ml
<b>Target</b>	YTHDF3
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
<b>Uniprot ID</b>	<a href="#">Q7Z739</a>
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.



Human YTHDF3(YTH domain-containing famil...