


Product Datasheet

Human LACTB2(Endoribonuclease LACTB2) ELISA Kit (orb1715175)

Description	Human LACTB2(Endoribonuclease LACTB2) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	LACTB2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q53H82
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
LACTB2(Endoribonuclease
LACTB2) EL...