


## Product Datasheet

# Porcine ADP/Acrp30(Adiponectin) ELISA Kit (orb1715196)

<b>Description</b>	Porcine ADP/Acrp30(Adiponectin) ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	3.125-200pg/ml
<b>Target</b>	ADP/Acrp30
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<a href="#">Q6PP07</a>
<b>Sensitivity</b>	1.875pg/ml
<b>Expiration Date</b>	Please enquire.



Porcine  
ADP/Acrp30(Adiponectin)  
ELISA Ki...