


Product Datasheet

Human EHMT1(Histone-lysine N-methyltransferase EHMT1) ELISA Kit (orb1715197)

Description	Human EHMT1(Histone-lysine N-methyltransferase EHMT1) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/ml
Target	EHMT1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9H9B1
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
EHMT1(Histone-
lysine N-
methyltrans...