

## Product Datasheet

# Human Perforin Monoclonal Antibody, FITC Conjugated Clone B-D48 (orb1715599)

<b>Description</b>	Human Perforin Monoclonal Antibody, FITC Conjugated Clone B-D48
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	FC
<b>Immunogen</b>	Natural human Perforin
<b>Target</b>	Recognizes human Perforin, a 70 kDa protein
<b>Form/Appearance</b>	Phosphate-buffered saline with 5% rAlbumin and 0.09% sodium azide
<b>Storage</b>	Stable at +2-8°C for 12 months. For longer storage freeze aliquots
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Source</b>	Myeloma X63/AG.8653 x Balb/c spleen cells
<b>Dilution Range</b>	Use 10 µl to label 10 <sup>6</sup> cells or 100 µl of whole blood
<b>Expiration Date</b>	12 months from date of receipt.