

Product Datasheet

Human IL-34 Detection Antibody, Biotin Conjugated Clone B-R34 (orb1715982)

Description	Human IL-34 Detection Antibody, Biotin Conjugated Clone B-R34
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Note	For research use only
Application notes	This biotinylated Antibody can be used as Detection in a human sandwich Immunoassay. The suggested coating concentration range should be optimised by each laboratory for each application.
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.