

## Product Datasheet

**TIA-1 Antibody / T-Cell Intracellular Antigen-1  
(orb1716151)**

**Description**

The T cell intracellular antigen 1 (TIA

associated protein also d...

**Species/Host**

Mouse

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested Applications**

IHC-P

**Immunogen**

Recombinant full-length human protein was used as the immunogen for the T-Cell Intracellular Antigen-1 antibody.

**Storage**

Aliquot the T-Cell Intracellular Antigen-1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

**Note**

For research use only

**Application notes**

Optimal dilution of the T-Cell Intracellular Antigen-1 antibody should be determined by the researcher.

**Formula**

1 mg/ml in 1X PBS; BSA free, sodium azide free

**Isotype**

Mouse IgG1, kappa

**Clonality**

Recombinant

**Uniprot ID**

**P31483**

**Hazard Information**

This T-Cell Intracellular Antigen-1 antibody is available for research use only.

**Dilution Range**

Immunohistochemistry (FFPE): 1-2ug/ml

**Expiration Date**

12 months from date of receipt.