
Product Datasheet

**TIA-1 Antibody / T-Cell Intracellular Antigen-1
(orb1716153)**

Description

The T cell intracellular antigen 1 (TIA

associated protein also d...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

IHC-P

Immunogen

Recombinant full-length human protein was used as the immunogen for the T-Cell Intracellular Antigen-1 antibody.

Storage

Aliquot the T-Cell Intracellular Antigen-1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Optimal dilution of the T-Cell Intracellular Antigen-1 antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Recombinant

Uniprot ID

[P31483](#)

Hazard Information

This T-Cell Intracellular Antigen-1 antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml

Expiration Date

12 months from date of receipt.