



Product Datasheet

MGMT Antibody (orb1716200)


Description	Cancer chemotherapeutic alkylating agents (e.g. BCNU,) act by inducing formation of lethal cross lin...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Purified His-tagged MGMT protein was used as the immunogen for the MGMT antibody.
Storage	Aliquot the MGMT antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the MGMT antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	P16455
Hazard Information	This MGMT antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human lymph node ti...



IHC staining of FFPE human stomach tissu...



IHC staining of FFPE human testis tissue...